

BASELINE ISSUES Causes Mobile phase degrade Degraded mobile phase can cause noisy baseline. **Noisy Baseline** - Prepare fresh mobile phase. Flush system and column with fresh mobile phase. NOISE Inappropriate grade mobile phase - Use HPLC/UHPLC or LC-MS grade solvent. Air in mobile phase flow line - Air in flow line can cause noisy baseline or random - Degas mobile phase. Purged pump using degassed mobile phase. - Use mixer with appropriate volume. Consult your manufacturer for more details on therecommended mixer volume

	Causes	Solutions
Baseline Drift DRIFT	Mobile phase absorption	Base line will drift when using UV absorbing solvent as mobile phase, particularly in gradient elution. Ensure detection wavelength is above UV cut-off of solvent. Use alternative solvent that does not absorbs UV at detection wavelength.
	Contaminant in mobile phase	- Prepare fresh mobile phase

pikes	Causes	Solutions
· 	Air in mobile phase flow line	- Air in flow line can cause noisy baseline or random spikes. - Degas mobile phase. - Purged pump using degassed mobile phase.
- 	Pump Issue	- Contact service engineer for assistance.

	ı	PRESSURE ISSUES
	Causes	Solutions
	Buffer precipitation	- Ensure compatibility of buffer-organic mixture to prevent buffer precipitation In the event of accidental buffer precipitation, flush system and column with low organic content mobile phase to dissolve the salt.
High Pressure	Sample precipitation	- Ensure sample is soluble in method's initial mobile phase composition to prevent precipitation of sample. - Contact service engineer if system clogged. - Flush column at half the flow rate with stronger solvent that can dissolve the sample if column is clogged.
₽ ^	Contaminated or clogged column	- Check if system is clogged by disconnecting the column. Contact service engineer if system is clogged If using guard column, check if guard column is clogged. Replace the guard column if clogged Check if main column is clogged or contaminated. Flush column at half the flow rate with stronger eluting solvent to remove contaminants. If pressure does not return to normal, reverse flush the column (if allowed) to remove particulate clogging the column inlet. If pressure still does not return to normal, replace the column Install an in-line filter in front of analytical column to prevent clogging of column due to particulates Use a guard column to prevent contamination of column by strong binding compounds.

Colution

Carra

flow line

check vale.

Faulty pump seal or

LOW OI	Causes	SOLUTIONS
No Pressure	Insufficient mobile phase	Replace mobile phase and ensure it is sufficient for analysis. Purge pump and system with mobile phase to remove air in flow line. Flush column with mobile phase. If column performance is affected, replace column.
	Mobile phase leak	- Check for leak in mobile phase flow line Ensure fitting is sufficiently tight to prevent leakage Use fitting with appropriate pressure rating to prevent leakage.
Unstable	Causes	Solutions
Pressure	Dubbles in mobile whose	B 13 1

Purge pump to remove bubbles.

- Contact service engineer for assistance

LC TROUBLESHOOTING TIPS



BASIC STEPS Follow these three steps to isolate where the problems is. Check the obvious explanations first and change only one thing at

Check the Basics:

- Power supply
 Analytical conditions
 Electrical connections
- Mobile phase volume
 Mobile phase leakage

Identify the Cause:

- Define the problem clearly; for example, "Peak tailing occurs after changing
- Review sample preparation procedure and maintenance records to identify trends in the data or problem indicators, such as pressure increasing over time.
- Use a logical sequence of steps to isolate possible causes.

Document Everything:

- Document all troubleshooting steps and results; this may help you identify and solve the next problem faster.
- Always inject a test mix and compare to previous data to ensure restored performance.

Still having problems? Let us help. Contact us at: consumablesap@shimadzu.com.sg



Causes - Check for leak in mobile phase flow line. - Ensure fitting is sufficiently tight to prevent Mobile Phase Leak leakage - Use fitting with appropriate pressure rating to prevent leakage. Wrong Detector Ensure correct data acquisition time setting. Ensure correct wavelength setting for UV or No Peak florescence detector. - Ensure correct mass setting for MS detector. - Ensure correct vial number setting. - Ensure sufficient sample volume in vial. No Injection Ensure injector needle not clogged. Sample Concentration Increase sample concentration or injection - Increase elution strength of mobile phase to elute compound. Sample Not Eluted

Solutions

activity.

original performance.

- Use endcapped reversed phase columns to reduce silanol activity.

- Adjust pH of mobile phase to reduce silanol

- Flush column in reversed direction if allowed (check with column manufacturer for details).
- Replace column if flushing cannot recover

- Adjust mobile phase pH (do not exceed column pH limit).

Causes

group

Sample interacting with silica silanol

Inappropriate mobile phase pH Clogged column




Shop for

consumables

now at

ShopShimadzu.com

Column damaged - Replace column or worn out Causes Solutions Peak Fronting Sample overload - Decrease sample concentration or reduce niection volume Inappropriate diluent - Use mobile phase as diluent. Reduce injection volume Column damaged Flush column in reversed direction if allowed (check with column manufacturer for details). Replace column if flushing cannot recover original performance.



	Causes	Solutions
Change in Response	Sample issue	- Check sample concentration Check sample preparation procedure Check sample decomposition/shelf life.
	Mobile phase leak	- Check for leakage in flow line after injection port. - Ensure fitting is sufficiently tight to preventleakage. - Use fitting with appropriate pressure rating to prevent leakage.
	Detector lamp worn out	- Check that detector lamp does not exceed recommended operating hours. - Replace detector lamp if exceed

recommended operating hours.



Causes

PEAK ISSUES


	Causes	Solutions
	Inappropriate diluent	- Use mobile phase as diluent.
	Contaminated Column	- Flush column with stronger eluting solvent to remove contaminants. - If problem persist, replace column.
=	Column Partially Clogged	- Reverse flush column with lower flowrate if allowed (check with column manufacturer for details).
	Column Damaged	- Replace column.

Solutions

	Mobile phase leak	Mobile phase leak will cause retention tim to increase. Check for leak in mobile phase flow line. Ensure fitting is sufficiently tight to preven leakage Use fitting with appropriate pressure rating to prevent leakage.
	Variation in column temperature	- Use a column oven to ensure consistent column temperature.
	Inadequate column equilibration	- Increase equilibration time.
arying tention Time	Evaporation of mobile phase	- Evaporation of mobile phase can cause change in the mobile phase composition. - Use solvent bottle cap to reduce solvent evaporation.
\mathcal{M}	Sample overload	- Decrease sample concentration or reduce injection volume.
	Column stationary phase collapse/ dewetting	- Phase collapse/dewetting will cause retention time to decrease Rewet column using high organic content mobile phase Replace column if rewetting fails If analysis require use of low organic content mobile phase, ensure the column used is compatible with 100% aqueous condition.
	Degradation of stationary phase	- Degradation of stationary phase will cause retention time to decrease. - Replace damaged column. - Ensure mobile phase pH within recommended pH limit.
	Pump issue	- Check flow rate of each pump. - If flow rate is lower than setting, contact service engineer for assistance.
	Causes	Solutions

_	Causes	Solutions
Poor esolution	Non-selective stationary phase	- Choose appropriate stationary phase and column dimension.
	Poor efficiency	Increase column length.Decrease column particle size.
	Sample overload	- Optimise sample concentration or sample injection volume.
	Incorrect analysis conditions	- Check and optimise gradient program, flow rates and column temperature etc.

Carryover/ Ghost Peak	Causes	Solutions
	Previous analysis ended too soon	- Increase analysis run time Increase elution strength of mobile phase to elute stronger retained compounds.
	Contaminated mobile phase	 Perform analysis without injection. If peak is observed, contaminant originated from mobile phase. Use appropriate grade solvent (Organic: HPLC/ UHPLC/LC-MS grade, aqueous: ultrapure water). Prepare fresh mobile phase daily.
		- Use Ghost Trap (For non-MS detector) or GLC Suction filter 2. (contact Shimadzu for more information.)
	Contaminated diluent	 Perform diluent blank analysis. If peak is present, contaminant originated from diluent. Use appropriate grade diluent (Organic: HPLC/ UHPLC/LC-MS grade, aqueous: ultrapure water).
	Contaminated injection port/	Perform needle and injection port rinse. Use appropriate solvent for needle and injection.







Install **Inline filter** between auto-sampler and column to remove pump seal debris and sample particles. This will prevent column back pressure problems.

Have pressure issues?

particulates.



Want to eliminate Ghost Peaks?

Use **Ghost Trap DS** to effectively adsorb impurities in the mobile phase in order to reduce the time required for method development and impurity analysis.

Use **GLC Suction Filter 2** to remove contaminants in mobile phase to minimizes chromatographic interferences. It is suitable for LC/MS analysis due to its low bleed.

