

AOAC 2023 TO67

Nozomi Maeshima<sup>1</sup>, Dominika Gruszecka<sup>2</sup>, Manami Kobayashi<sup>1</sup> 1 Shimadzu Corporation, Kawasaki, Kanagawa, Japan; 2 Shimadzu Scientific Instruments Inc., Columbia, Maryland



## Introduction

Per- and Polyfluoroalkyl Substances (PFAS) is the collective name for a chemical group of organic fluorinated compounds, perfluorooctane sulfonate (PFOS) and perfluorooctanoic acid (PFOA) are representative compounds of PFAS. They have been used water repellents, surface treatment agents, fire extinguishers, and coatings. PFAS are persistent and bioaccumulative in the environment because of their stable structure and known that they are present in a wide range of environmental water and wildlife. Due to concerns about human exposure through diet, studies on the status of food contamination by PFAS are being conducted in various countries. We have examined a quantitative analysis method for forty PFAS compounds in foods with two SPE cartridges.

### Methods and Results

#### Material, sample, and equipment

Standard compounds were purchased from Wellington Laboratories. Carrot was purchased from a grocery store and homogenized using a freeze grinder. Quantification was performed using fully polypropylene low-binding vials TORAST™-H Bio Vial (Shimadzu GLC Ltd.) with a triple

quadrupole mass spectrometer LCMS-8050 equipped with Nexera™ X3 UHPLC (Shimadzu). The system configuration is shown below. To prevent contamination from an equipment, a delay column was added between a mixer and an autosampler.

#### Nexera X3 system

Column : Shim-pack Scepter™ C18-120 (100 mm x 2.0 mm I.D., 1.9 µm)
Delay column : Shim-pack GIST C18 (3.0 mm x 50 mm I.D., 5 µm, HSS)

Mobile phase A: Acetonitrile/water=5:95(v/v) containing 2 mmol/L Ammonium acetate

Mobile phase B: Acetonitrile

Rinse : Methanol/water=50:50(v/v)

Flow rate : 0.4 mL/min

Time program : B conc. 20% (0 min)  $\rightarrow$  100% (10-12 min)  $\rightarrow$  20% (12.01-15 min)

The flow was introduced into the mass spectrometer between 0.1 to 9.6 min using a

flow switching valve.

Column temp. : 40 °C Injection vol. : 1  $\mu$ L

#### LCMS-8050

Ionization: ESI, Negative modeDL temp.: 200 °CInterface temp.: 200 °CHeat block temp.: 300 °CNebulizer gas: 3 L/minHeating gas: 10 L/minDrying gas: 10 L/minProbe position: +2 mm



#### Extraction

Extraction was performed using a pre-processing method, taking reference from the QuEChERS method. It is a simple protocol that does not require glassware, and

centrifugation is required only once. The procedure is shown in figure 1.

#### Development of purification process

Initially, purification was performed refer to 2nd draft method 1633 of EPA. However, due to significant losses in purification step with the SPE cartridge, detailed investigation was conducted focusing on the washing and elution steps. Two types of weak anion-exchanged SPE cartridges were evaluated: InertSep MA-2 (GL Sciences Inc.) or EVOLUTE® EXPRESS WAX (Biotage). Eleven aqueous solutions for washing and elution were prepared with methanol concentrations ranging from 0 to 100% in 10% increments. Three types of eluents were prepared: one containing formic acid (formic acid/methanol solution=1:1000(v/v/)), one containing ammonium (28% ammonia solution/methanol solution=1:100(v/v/v)), and one containing nothing, and one set of eluent was

sequentially supplied to each SPE cartridge. After loading the carrot extract containing PFAS into the SPE cartridge, elution was performed starting from a 5 mL portion with 0% methanol ratio, then gradually increasing (figure 2). Biotage® PRESSURE+ 48 was used for this process.

The compounds present in each fraction were quantified (figure 3). Since there was a compound that eluted up to 53.5% (3:3 FTCA) at methanol ratio of 50%, formic acid/methanol/water=1:400:600(v/v/v) was chosen as the washing solution. Only a small number of compounds required methanol ratio of 100% for elution, so 28% ammonia solution/methanol/water=1:90/10(v/v/v) was used as the eluent.

#### **Evaluation of concentration**

The feasibility of concentrating the eluate using nitrogen gas blow-down was investigated. After adding standard compounds to 5 mL of purified solution of carrot, it was concentrated to less than 1 mL and filled up to 1 mL with methanol. Upon quantification with LC-MS, there were loss ranging from 67.6 to 98.7% for eight compounds,

PFOSA, NMeFOSA, NEtFOSA, NMeFOSE, NEtFOSE, 3:3 FTCA, 5:3 FTCA, and 7:3 FTCA. The losses for the other compounds were 27.8% or less. Therefore, the concentration step using nitrogen gas blow-down was avoided in this pre-processing method.






Figure 1. The extraction process




Figure 2. The purification process

Total of 12 fractions, including the through fraction,
were collected and individually analyzed to quantify PFAS
concentration present in each fraction.






Figure 3. Amount of elution at each methanol concentration
The black line in the figure is drawn between methanol ratios of 50% and 60%.






Figure 4. The purification process

#### Recovery rate test

To confirm the matrix effects, the eluate obtained carrot extract purification was analyzed after adding standard compounds (1 ng of PFOS in 5 mL of eluate). Additionally, to evaluate purification efficiency, standard compounds were added to the solution obtained from carrot extraction (1 ng of PFOS in 10 mL extraction 5-fold diluted with water), and the purified eluent was analyzed. Furthermore, to quantify the total loss in extraction step and purification step, ground carrot with added standard compounds (5 ng of PFOS in 10 g of carrot) were extracted, and the eluent was analyzed with LC-MS.

The recovery rates were calculated by comparing the peak areas of standard solutions, without correction by surrogate compounds. The results of recovery rate and Final concentration in vial before addition of formic acid are shown in table 1, and MS chromatograms are shown in figure 5. We obtained good recovery rate of within 70% to 120% for both cartridges. When injecting 5  $\mu$ L, 6:2FTS exhibited higher concentration than that of theoretically expected due to matrix effect. However, reducing the injection volume suppressed matrix effect, suggesting a good recovery rates for all compounds.






Figure 5. MS chromatograms of PFAS 40 compounds



Table 1. Recovery rate of PFAS using InertSep MA-2 or WAX

Compound		Conc. in carrot (ng/g)	Final conc. (ng/mL)	Recovery rate (%)					
	Retention time			InertSep MA-2			WAX		
	(min)			After purification	After extraction	Before extraction	After purification	After extraction	Before extraction
PFBA	1.38	2.0	0.8	93.3	105.1	100.2	89.7	106.5	106.3
PFPeA	2.55	1.0	0.4	90.5	100.9	96.2	86.8	100.0	100.7
PFHxA	3.33	0.5	0.2	97.7	102.7	100.8	97.7	101.3	109.4
PFHpA	3.92	0.5	0.2	101.0	107.8	96.9	93.4	105.0	106.4
PFOA	4.40	0.5	0.2	86.0	99.7	88.6	84.2	103.1	88.5
PFNA	4.84	0.5	0.2	93.7	103.3	102.5	97.6	108.7	104.6
PFDA	5.25	0.5	0.2	97.5	94.4	98.4	87.2	95.6	99.8
PFUnA	5.65	0.5	0.2	87.4	95.0	92.1	81.5	91.3	92.2
PFDoA	6.05	0.5	0.2	95.9	101.1	102.1	89.1	96.3	99.0
PFTrDA	6.43	0.5	0.2	82.3	83.5	89.4	79.2	82.5	87.3
PFTeDA	6.81	0.5	0.2	88.2	89.8	91.9	85.9	95.7	94.9
PFBS	3.39	0.5	0.2	100.5	94.0	103.7	104.5	105.6	105.9
PFPeS	4.04	0.5	0.2	98.1	98.1	104.7	96.3	96.9	91.7
PFHxS	4.56	0.5	0.2	100.7	88.7	92.1	85.4	98.5	86.4
PFHpS	5.03	0.5	0.2	92.3	104.6	100.9	97.1	105.9	102.8
PFOS	5.45	0.5	0.2	89.8	86.2	88.9	75.9	80.4	98.8
PFNS	5.87	0.5	0.2	95.6	97.7	95.6	94.9	96.2	98.7
PFDS	6.27	0.5	0.2	87.4	101.0	95.0	83.1	92.1	97.6
PFDoS	7.04	0.5	0.2	86.6	94.1	92.8	92.9	91.6	94.8
4:2FTS	3.11	2.0	0.8	78.0	93.5	92.9	81.6	99.6	99.1
6:2FTS	4.20	2.0	0.8	86.5	103.3	93.9	93.8	101.9	92.3
8:2FTS	5.05	2.0	0.8	85.7	92.9	91.9	84.3	90.0	97.7
PFOSA	7.23	0.5	0.2	95.0	92.5	98.5	91.3	93.5	100.0
NMeFOSA	8.76	0.5	0.2	91.2	99.8	92.1	86.3	94.6	104.6
NEtFOSA	9.13	0.5	0.2	92.0	102.8	100.3	90.9	92.1	106.3
NMeFOSAA	5.27	0.5	0.2	84.6	109.3	101.8	97.9	93.8	98.4
NEtFOSAA	5.45	0.5	0.2	104.5	111.4	79.0	90.6	99.9	114.0
NMeFOSE	8.60	5.0	2.0	91.4	94.2	95.9	85.2	88.7	99.3
NEtFOSE	8.98	5.0	2.0	93.8	94.8	99.0	87.7	89.2	98.9
HFPO-DA	3.56	2.0	0.8	87.8	99.3	98.0	85.1	99.6	103.5
ADONA	4.10	2.0	0.8	90.9	100.1	98.9	93.2	99.9	95.5
9Cl-PF3ONS	5.76	2.0	0.8	95.0	100.2	97.8	95.7	100.1	101.1
11Cl-PF3OUdS	6.57	2.0	0.8	92.8	97.4	96.0	93.3	96.1	100.1
3:3 FTCA	1.87	2.5	1.0	84.7	86.4	84.6	85.2	80.8	89.9
5:3 FTCA	3.63	12.5	5.0	91.2	100.1	93.6	87.9	98.0	98.7
7:3 FTCA	4.76	12.5	5.0	85.0	97.3	96.8	85.9	98.0	99.7
PFEESA	3.70	1.0	0.4	91.5	97.7	95.1	91.2	96.8	101.0
PFMPA	1.79	1.0	0.4	88.2	100.5	96.4	86.2	102.2	101.2
PFMBA	2.84	1.0	0.4	86.3	98.6	92.5	88.9	97.5	99.0
NFDHA	3.26	1.0	0.4	82.7	93.6	93.8	85.2	100.3	92.9
<70%				0	0	0	0	0	0
70-120				40	40	40	40	40	40
>120%				0	0	0	0	0	0



#### Discussion

Purification conditions using SPE cartridges were studied and the optimized eluents was selected. Purification was performed using the optimized process and good recoveries were obtained. Two SPE cartridges of InertSep NA-2 and WAX were compared, but no significant difference was observed.

## Conclusions

- An LC-MS method for forty PFAS within fifteen minutes analysis were created.
- The development of the pre-processing step, and a recovery test were conducted under low PFAS concentration conditions, resulting in favorable results.

## References

Analysis of Per- and Polyfluoroalkyl Substances (PFAS) in Aqueous, Solid, Biosolids, and Tissue Samples by LC-MS/MS, 2nd Dragt Method 1633, EPA (June 2022).

The products and applications in this presentation are intended for Research Use Only (RUO). Not for use in diagnostic procedures.



Shimadzu Corporation www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

First Edition: October, 2023