

# DIOXIN 2013 P-0001

Teaseob Shin<sup>1</sup>, Hyunook Kim<sup>1</sup>, Youngmin Hong<sup>2</sup>, Insun Lee<sup>2</sup> <sup>1</sup>The University of Seoul, <sup>2</sup>Dong-il Shimadzu Corporation



# 1. Introduction

Triclosan and triclocarban are antimicrobials currently added to a large number of consumer's products including household soaps, deodorants, toothpastes, mouth washes, and cleaning agents. Such release of triclosan and triclocarban has the potential to cause a number of environmental and human health problems. Triclosan is regulated all by the U.S. Food and Drug Administration, the Environmental Protection Agency, and the European Union. By wastewater treatment processes, a portion of triclosan can be degraded and partly the remaining adsorbs to sewage sludge or biosolids. In general, analysis for quantification of triclosan and triclocarban in water and solid samples has been used LC/MS/MS through the pretreatment and preconcentration. In this study, an on-line column switching pre-concentration method in which a surface modified polymer resin is used has been applied to extract triclosan, triclosan-methyl, and triclocarban in water samples. The pre-extracted chemicals were quantified using LC/MS/MS with large volume injection. The large volume injection enhances not only the LC/MS/MS throughput, but also the trapping efficiency and recovery of the analysis.

# 2. Experimental

## 2-1. Physico-chemical properties



### 2-2. Analytical condition

1) HPLC : Shimadzu HPLC system

1	/ TIFEC . SHITIAUZU TIFEC System					
	Pre-concentration	: Shim-pack MAYI-ODS (G) (4.6 mml.D. × 10 mmL., 50 µm)				
	Separation column	: ACE3 C18-AR (3.0 mml.D. × 100 mmL., 3 μm)				
	Mobile phase A / B	: 5 mM ammonium formate / methanol				
	Mobile phase C	: 10% Methanol in water				
	Gradient program	: 70%B (5.5 min) -100%B (8-10 min) -70%B (10-20 min)				
	Flow rate A+B	: 0.5 mL/min				
	Flow rate C	: 1.0 mL/min				
	Injection volume	: 1500 μL				
	Column temperature	: 40°C				

2) MS : Shimadzu LCMS-8040 Triple quadrupolemass spectrometer

: 3 L/min
: 15 L/min
: 250°C
: 400°C
: ESI negative ion mode
: MRM mode (multiple reaction monitoring)
Pause time 3 ms, Dwell time 20 ms



	lon	Precursor ion	Product ions (CE, V)	
Triclosan	[M-H] <sup>-</sup>	287	35.2 (9)	142.1 (35)
Triclocarban	[M-H] <sup>-</sup>	313	160.1 (14)	126.2(20)





\* Sample preparation – sample 1.7 mL + methanol 0.3 mL

Fig. 1 Diagram of system flow channel

# 2-3. Schematic diagram of sewage treatment plant under study



\* Sampling points of the sewage treatment process

# 3. Results & Discussion

### 3-1. Filtration test

(n =10, ng/L)

		Tricle	osan	Triclocarbar			
Material	PVDF	PTFE	Cellulose Acetate	PVDF	PTFE	Cellulose Acetate	
Pore size	0.2 µm	0.2 µm	0.2 µm	0.2 µm	0.2 µm	0.2 µm	
Before filtration		997.1			999.6		
After filtration	591.7	887.6	493.2	616.8	925.2	514.0	

### 3-2. Linearities of calibration curves (n=7)





# 3-3. LCMS Chromatogram



### 3-4. Method validation

Table 1 Recovery test							
			,			(n = 7)	
		DW	Influent	Effluent	MDL(ng/L)	C.V.(%)	
Telelocop	0.1 ng/mL	98%	91%	92%	2.0	9.1	
Inclosen	1 ng/mL	99%	94%	91%	2.8		
Triclocarban	0.1 ng/mL	98%	91%	91%	2 5	0 1	
ITICIOCAIDAN	1 ng/mL	99%	89%	91%	2.0	0.1	

### 3-5. Sample analysis

Table 2 Sewage treatment plant water sample collected from each process

ng/L	Influent	Primary clarifier	Anaerobic	Anoxic	Aerobic	Secondary clarifier	Effluent
Triclosan	75	109	63	77	48	31	28
Triclocarban	72	75	48	39	47	45	38

# 4. Conclusion

The developed method using a methylcellulose immobilized pre-concentration column switching system could be applied successfully to quantitate residual triclosan and triclocarban without pretreatment in wastewater samples. The preliminary finding suggests that pre-concentration method is potentially useful for on-line pre-concentration method for analysis of trace level triclosan and triclocarban in wastewater samples as well as other application.

First Edition: October, 2013



Shimadzu Corporation www.shimadzu.com/an/

#### For Research Use Only. Not for use in diagnostic procedures.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.