

ASMS 2013 ThP 115

Shannon L. Cook¹, Hideki Yamamoto², Tairo Ogura², Yusuke Osaka², Ichiro Hirano²
¹Shimadzu Scientific Instruments. 7102 Riverwood Drive Columbia, MD 21046, USA 2Shimadzu Corporation.
1, Nishinokyo-Kuwabaracho Nakagyo-ku, Kyoto

604–8511, Japan



1. Introduction

Recently, the promotion of proteomic analysis has resulted in the identification of many proteins. Analyzing the proteins function, such as its various signaling pathways, is a clue to understanding disease. Historically, quadrupole time-of-flight mass spectrometry has been preferred over more quantitative, but less sensitive, triple quadrupole based mass spectrometric methodologies. Our aim is to develop a triple quadrupole based system to quantitatively analyze peptide fragments through a shotgun proteomics approach. We developed both a conventional and nano-flow quantitative system for proteomic analyses utilizing liquid chromatography (LC) coupled to a triple quadrupole mass spectrometer (triple Q MS). We also demonstrate the sensitivity enhancement of the triple Q MS as a function of the MS front-end LC at analytical and nano-flow rates (Fig. 1).



Fig. 1 Nano-flow LC and LCMS-8040 triple quadrupole mass spectrometer

2. Methods and Materials 2-1. Methods

Samples were measured by ESI-MS coupled to either an UHPLC system or a nano-flow HPLC system.

Analysis by UHPLC coupled to triple Q MS (conventional LC-MS)

UHPLC conditions (Nexera system)

Column	: Shim-pack XR-ODS III 2 mm I.D. $ imes$ 50 mm L., 1.6 μ m						
Mobile phase A	: 0.5% acetic acid in aq.						
В	: acetic acid / Water / Acetonitrile (0.5/19.5/80)						
Flow rate	: 0.2 mL/min						
Time program	: B conc.5%(0 min) - 60%(10 min) - 100%(10.01-20 min) - 5%(20.01-30 min)(Fig. 2)						
Injection vol.	: 1 µL						
Column temperature	: 40°C						

MS conditions (LCMS-8040)

Ionization	: ESI
Events	: Positive MRM mode
Applied Voltage	: 4.5 kV
Desorption Liquid temp.	: 300°C
Heat Block temp	: 500°C
Nebulizer Gas flow rate	: 3 L/min
Drying Gas flow rate	: 10 L/min



Fig. 2 The time program of conc.B gradient curve of the UHPLC system



Analysis by nano-flow HPLC coupled to triple Q MS (nano-flow LC-MS)

nanoLC conditions (Pr	rominence system)					
Column	: column with nano spray tip 100 μ m I.D. $ imes$ 150 mm L., 3 μ m (Reprosil ^R particles)					
Mobile phase A	: 0.5% acetic acid in aq.					
В	: acetic acid / Water / Acetonitrile (0.5/19.5/80)					
Flow rate	: 500 nL/min					
Time program	: B conc.5%(0 min) - 60%(30 min) - 100%(31-40 min) - 5%(45-65 min) (Fig. 3)					
Injection vol.	: 1 µL					
Column temperature	: Room Temp.					

MS conditions (LCMS-8040)

Ionization	: nano-ESI (sprayed by nano spray tip)
Events	: Positive MRM mode
Applied Voltage	: 2.5 kV
Desorption Liquid temp.	: 250°C
Heat Block temp.	: 200°C
Nebulizer Gas flow	: None
Drying Gas flow	: None



Fig. 3 The time program of conc.B gradient curve of the nanoLC system

The nano-flow HPLC system consists of binary pumps and autosampler. The flow channel is illustrated below (Fig. 4).



Fig. 4 The flow channel configuration of the nano-flow HPLC system

2-2. Materials

Digested BSA (Bovine Serum Albumin) was analyzed to optimize the conditions of the LC and the triple Q MS parameters.



Digestion of BSA

1 mg/mL of BSA in 50 mM Ammonium Bicarbonate (ABC) and 8 M Urea solution
↓ Add 10 mM of dithiothreitol, 37°C 30 min
↓ Add 50 mM of 2-iodoacetamide, 37°C 30 min
↓ Add 10 µL of LysC (1 µg/µL), 37°C 4h
↓ Diluted by 3 mL of 50 mM ABC
↓ Add 20 µg of Trypsin, 37°C overnight
Quenched by 10% trifluoroacetic acid (TFA) in aq.

Desalination of peptides

Add 1 mL of Buffer B to Stage Tip ↓ Centrifuge 500 g 3 min Add 1 mL of Buffer A to Stage Tip ↓ Centrifuge 500 g 3 min Add 1 mL of BSA Digestion to Stage Tip ↓ centrifuge 500 g 4 min Add 1 mL of Buffer A to Stage Tip ↓ centrifuge 500 g 4 min Add 250 µL of Buffer B to Stage Tip ↓ Centrifuge 300 g 2 min Preserve elution buffer at -20°C

3. Results

Most peptides in digested BSA were detected as protonated divalent or trivalent molecules([M+2H]²⁺ or [M+3H]³⁺) in ESI positive ion mode. The precursor ions were monitored by Q1 SIM measurement. Mass spectrometric parameters for MRM analysis, such as product ion and collision energy, were optimized by an automatic MS optimization procedure. The chromatographic and automatic MS optimization produced 107 MRM transitions for 11 peptides for both the nano-flow LC-MS method and the conventional LC-MS methods. Quantitative analyses were performed for peptides containing at least 7 amino acids (Table 1).



Calibration curves were generated from each MRM transition with linearity improving when analyzed by nano-flow LC-MS as opposed to conventional LC-MS. Concentration ranges of the conventional LC-MS method and the nano-flow LC-MS method were 1-50 pmol and 10-500 fmol. Additionally, the LOD and LOQ of each transition analyzed by the nano-flow LC-MS method is several hundred fold lower than the conventional LC-MS method (Table 1, Fig. 6). In summary, we improved the sensitivity of a triple Q MS by reducing the LC flow rate to nL/min. A few femtomoles of injected peptides could be analyzed by nano-flow LC-MS.



Fig. 6 Representative calibration curve [DDSPDLPK>b2(443.70>230.90)]



4. Conclusions

- Nano-flow LC enhanced the sensitivity of a triple Q MS, as a function of the MS front-end LC, as opposed to conventional LC while improving the linearity of calibration curves.
- We developed and evaluated a novel quantitative system for proteomics with the results suggesting that this system could be applicable as a tool for shotgun proteomics.

Table 1 MRM transitions of BSA Digests and the LOD, LOQ, and relation coefficient (R²) of each transition for both the nano-flow LC-MS method and the conventional LC-MS method.

		conventional LC-MS			nano-flow LC-MS		
MRMtransition	m/z	LOD(pmol)	LOQ(pmol)	coefficient(R2)	LOD(fmol)	LOQ(fmol)	coefficient(R2)
CCDKPLLEK>y2++-H2O	646.30>129.30	0.41	1.24	0.996	0.58	1.74	1.000
CCDKPLLEK>y1	646.30>147.00	0.36	1.08	0.999	1.67	5.02	1.000
CCDKPLLEK>y10++-H2O	646.30>637.35	0.08	0.24	0.995	0.27	0.80	0.999
CCDKPLLEK>y3+++	646.30>130.20	0.29	0.86	0.996	3.92	11.75	0.999
CCDKPLLEK>b4+++-H2O	646.30>183.30	0.30	0.91	0.997	2.25	6.75	0.999
CCDKPLLEK>y2-H2O	646.30>258.10	0.19	0.56	0.997	1.94	5.83	0.995
CTESLVNR>b4+++-H2O	569.75>178.00	0.03	0.10	0.998	0.41	1.24	1.000
CTESLVNR>y9++-H2O	569.75>560.95	0.03	0.09	0.998	1.20	3.61	0.999
CTESLVNR>y7	569.75>818.45	0.02	0.06	0.997	0.32	0.97	0.999
CTESLVNR>y1	569.75>175.20	0.03	0.09	0.996	1.91	5.72	0.998
CTESLVNR>b2	569.75>320.85	0.06	0.17	0.994	0.37	1.12	1.000
CTESLVNR>YZ	569.75>289.05	0.09	0.28	0.999	1.58	4.73	1.000
CTECLUMP	569.75>588.15	0.06	0.17	0.996	0.29	0.88	0.998
CTESLVINR>YD	569./5>/1/.15	0.05	0.15	0.996	0.58	1.74 E 43	0.998
CTECLARDA DE LA LA	560 75- 313 10	0.09	0.20	0.995	1.01	3.42	0.998
	443 705213.10	0.06	0.32	0.992	0.48	1.45	1.000
	443.70>230.30	0.00	0.17	0.000	0.40	0.88	1.000
	443.70>656.35	0.00	0.16	0.998	0.07	0.00	1.000
	443.70>030.33	0.05	0.15	0.995	0.07	2 37	0.000
DSPDLPK>v5	443 70>569 20	0.15	0.45	0.999	0.75	0.82	1 000
DSPDLPK>v1	443 70>147 25	0.13	0.40	0.999	1.65	4 95	1 000
DSPDLPK>v6++-H2O	443 70>319 90	0.08	0.10	0.996	3.12	9.37	0.998
DSPDLPK>v5++	443 70>285 30	0.21	0.64	0.996	1 14	3 42	1 000
DSPDLPK>b3-H2O	443.70>300.20	0.12	0.36	0.993	2.22	6.66	0.999
DSPDLPK>y6+++-H2O	443.70>213.40	0.51	1.54	0.999	1.88	5.64	0.997
DSPDLPK>b4+++	443.70>139.10	0.30	0.89	0.998	3.33	9.99	1.000
LGEEHFK>y6	487.75>746.20	0.05	0.14	0.989	0.36	1.07	1.000
LGEEHFK>b2	487.75>229.05	0.05	0.15	0.988	0.27	0.82	0.999
LGEEHFK>b3+++	487.75>96.45	0.21	0.62	0.981	1.26	3.79	0.994
LGEEHFK>y2+++	487.75>98.45	0.86	2.57	0.965	2.46	7.37	0.997
LGEEHFK>y6+++	487.75>249.20	0.08	0.24	0.992	0.37	1.11	0.995
LGEEHFK>y6++	487.75>373.45	0.11	0.32	0.983	0.38	1.15	0.999
iACLLPK>y2	379.70>244.10	0.05	0.14	0.991	0.40	1.19	1.000
iACLLPK>b2	379.70>129.20	0.03	0.08	0.990	0.35	1.04	1.000
iACLLPK>y5++	379.70>315.70	0.04	0.12	0.990	0.37	1.11	1.000
iACLLPK>y5	379.70>630.35	0.04	0.13	0.990	0.38	1.13	1.000
iACLLPK>y1	379.70>147.20	0.04	0.12	0.992	0.45	1.34	1.000
ACLLPK>b2	379.70>289.10	0.04	0.12	0.990	0.39	1.16	1.000
ACLLPK>y3	379.70>357.15	0.05	0.14	0.991	0.32	0.95	1.000
ACLLPK>y4	379.70>470.40	0.05	0.14	0.993	0.39	1.16	1.000
ACLLPK>y6++	379.70>351.15	0.04	0.12	0.991	0.31	0.94	0.999
ACLLPK>b4	3/9./0>402.25	0.03	0.08	0.994	0.57	1./1	0.999
ACLLPK>D4++	3/9./0>201.25	0.03	0.10	0.992	0.68	2.03	0.999
VIDLIK>y5	395.25>5/7.35	0.03	0.09	0.991	0.37	1.10	1.000
VIDLIK>D2	395.25>212.90	0.04	0.11	0.991	0.38	1.14	1.000
	205 25> 249 05	0.02	0.05	0.989	0.50	1.51	0.000
VTDLTK>y2	205 25>246.05	0.03	0.14	0.989	0.40	1.20	0.999
VTDLTKSyf	205 25-676 25	0.04	0.02	0.991	0.04	0.53	0.999
VTDLTK>v3	395 25 361 10	0.05	0.08	0.991	0.17	0.52	0.999
VTDLTK v2-H20	395 25>230 10	0.05	0.75	0.990	0.25	1.04	0.998
FFVFVTK>h2	461 75>201 15	0.04	0.11	0.991	0.27	0.81	1 000
FFVFVTK>v6	461 75>722 40	0.03	0.10	0.992	0.41	1 24	1.000
EFVEVTK>v2	461 75>248 30	0.03	0.10	0.993	0.34	1.01	0.999
EFVEVTK>v5	461 75>575 30	0.05	0.16	0.994	0.29	0.88	0.999
EFVEVTK>y6++	461.75>361.60	0.06	0.17	0.991	0.25	0.74	0.999
EFVEVTK>y4	461.75>476.30	0.03	0.08	0.991	0.36	1.07	0.999
EFVEVTK>y3	461.75>346.95	0.04	0.11	0.989	0.26	0.77	0.998
EFVEVTK>y2-H2O	461.75>230.40	0.04	0.11	0.993	0.23	0.70	0.998
EFVEVTK>b2-H2O	461.75>183.00	0.04	0.12	0.996	0.22	0.65	0.997
EFVEVTK>y3+++-H2O	461.75>110.05	0.08	0.25	0.988	0.38	1.14	0.997
ACFAVEGPK>y1	554.25>147.25	0.05	0.16	0.996	1.06	3.19	0.996
ACFAVEGPK>y2	554.25>244.25	0.04	0.13	0.994	0.39	1.18	1.000
ACFAVEGPK>b3-H2O	554.25>342.95	0.04	0.12	0.991	0.36	1.09	1.000
ACFAVEGPK>y6	554.25>600.25	0.03	0.08	0.991	0.41	1.22	1.000
ACFAVEGPK>y3	554.25>301.00	0.04	0.12	0.993	0.50	1.49	0.999
ACFAVEGPK>y7	554.25>747.20	0.02	0.07	0.992	0.67	2.02	0.999
ACFAVEGPK>b2-H2O	554.25>182.95	0.05	0.16	0.992	0.49	1.47	0.999
ACFAVEGPK>y5	554.25>529.35	0.02	0.05	0.990	0.30	0.90	1.000
ACFAVEGPK>y4	554.25>430.05	0.03	0.10	0.989	0.35	1.06	0.999
ACFAVEGPK>y8++	554.25>454.25	0.04	0.12	0.991	0.39	1.17	1.000
ACFAVEGPK>y9++	554 25>490 05	0.02	0.06	0.988	0.34	1 02	0 999

First Edition: June, 2013



For Research Use Only. Not for use in diagnostic procedures.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

Shimadzu Corporation