

Shim-pack HPLC Column Selection Guide










Shim-pack Scepter

Wide usable pH range, Higher versatility for various applications, effective for method scouting

- Next-generation organosilica hybrid based columns with excellent chemical resistance to acids and alkalis at high temperatures.
- Strict column packing procedure that ensures robust and reproducible column performance.
- Available in a wide range of bonded chemistries and column dimensions.
- Shim-pack Scepter (metal-free) is an inert column with a PEEK-lined stainless steel body, ideal for analysis of metal-coordinating and ion-adsorbing compounds such as phosphate-containing and basic analytes.
- Shim-pack Scepter Claris features a column body with a newly-developed bioinert coating applied to the column body and stainless steel frit. Ideal for analysis of metal-coordinating and hydrophobically adsorbing compounds such as nucleic acids, proteins, and lipids.






 <div>Shim-pack Scepter C18-120</div>	<table><tr><td>Bonded Phase</td><td>Octadecyl</td></tr><tr><td>Carbon Loading</td><td>20%</td></tr><tr><td>Surface Area</td><td>360m²/g</td></tr><tr><td>Pore Size</td><td>120Å</td></tr><tr><td>Particle Size</td><td>1.9, 3, 5µm</td></tr><tr><td>pH Range</td><td>1~12</td></tr><tr><td>End-cap</td><td>Yes</td></tr><tr><td>USP Classification</td><td>L1</td></tr></table>	Bonded Phase	Octadecyl	Carbon Loading	20%	Surface Area	360m ² /g	Pore Size	120Å	Particle Size	1.9, 3, 5µm	pH Range	1~12	End-cap	Yes	USP Classification	L1	<div>General C18 Column:</div> <ul style="list-style-type: none">• Ideal for method development over a wide pH range of 1-12• Excellent chemical stability under a wide variety of acids, alkalis, and salts at temperatures up to 90 °C• High carbon load, higher separation power for complex samples	<div>Pharmaceutical Analysis</div> <div>Bioanalysis</div> <div>Food / Cosmetic / Environmental Analysis</div>	<ul style="list-style-type: none">• Macrolides antibiotics, -sartan drugs, -tidine drugs• Peptide therapeutics, oligonucleotides, ASO, siRNA• Natural products or TCM (e.g., alkaloids)• Generic drugs bioequivalence studies• Drugs and metabolites• TCM metabolites• Food preservatives, food dyes• Pesticides and veterinary drug residues• Cosmetic products (preservatives, steroids, etc)
Bonded Phase	Octadecyl																			
Carbon Loading	20%																			
Surface Area	360m ² /g																			
Pore Size	120Å																			
Particle Size	1.9, 3, 5µm																			
pH Range	1~12																			
End-cap	Yes																			
USP Classification	L1																			
 <div>Shim-pack Scepter C18-300</div>	<table><tr><td>Bonded Phase</td><td>Octadecyl</td></tr><tr><td>Carbon Loading</td><td>N.D.</td></tr><tr><td>Surface Area</td><td>N.D.</td></tr><tr><td>Pore Size</td><td>300Å</td></tr><tr><td>Particle Size</td><td>1.9, 3, 5µm</td></tr><tr><td>pH Range</td><td>1~12</td></tr><tr><td>End-cap</td><td>Yes</td></tr><tr><td>USP Classification</td><td>L1</td></tr></table>	Bonded Phase	Octadecyl	Carbon Loading	N.D.	Surface Area	N.D.	Pore Size	300Å	Particle Size	1.9, 3, 5µm	pH Range	1~12	End-cap	Yes	USP Classification	L1	<div>Wide pore C18 Column:</div> <ul style="list-style-type: none">• Ideal for method development over a wide pH range of 1-12• Excellent chemical stability under a wide variety of acids, alkalis, and salts at temperatures up to 90 °C• Large pore size suitable for analysis of medium to high molecular weight substances	<div>Pharmaceutical Analysis</div> <div>Bioanalysis</div>	<ul style="list-style-type: none">• Proteins, peptides, mAb, ADC• Oligonucleotides, mRNA• Generic drugs bioequivalence studies• Nucleic acid• Drugs and metabolites
Bonded Phase	Octadecyl																			
Carbon Loading	N.D.																			
Surface Area	N.D.																			
Pore Size	300Å																			
Particle Size	1.9, 3, 5µm																			
pH Range	1~12																			
End-cap	Yes																			
USP Classification	L1																			
 <div>Shim-pack Scepter HD-C18</div>	<table><tr><td>Bonded Phase</td><td>Octadecyl</td></tr><tr><td>Carbon Loading</td><td>25%</td></tr><tr><td>Surface Area</td><td>430m²/g</td></tr><tr><td>Pore Size</td><td>80Å</td></tr><tr><td>Particle Size</td><td>1.9, 3, 5µm</td></tr><tr><td>pH Range</td><td>1~12</td></tr><tr><td>End-cap</td><td>Yes</td></tr><tr><td>USP Classification</td><td>L1</td></tr></table>	Bonded Phase	Octadecyl	Carbon Loading	25%	Surface Area	430m ² /g	Pore Size	80Å	Particle Size	1.9, 3, 5µm	pH Range	1~12	End-cap	Yes	USP Classification	L1	<div>High Density C18 Column:</div> <ul style="list-style-type: none">• Superior stereoselectivity, effective in separation of structural isomers and analogs• Ultra-inert, use of adsorptive ion-pairing reagent (e.g. trifluoroacetic acid, triethylamine) should be avoided in mobile phase	<div>Pharmaceutical Analysis</div> <div>Bioanalysis</div> <div>Food / Cosmetic / Environmental Analysis</div>	<ul style="list-style-type: none">• Drug products, API, intermediates and impurities• Natural products or TCM, especially those containing alkaloids• Structural isomers and analogs• Food preservatives, food dyes• Cosmetic products (preservatives, steroids, etc)
Bonded Phase	Octadecyl																			
Carbon Loading	25%																			
Surface Area	430m ² /g																			
Pore Size	80Å																			
Particle Size	1.9, 3, 5µm																			
pH Range	1~12																			
End-cap	Yes																			
USP Classification	L1																			
 <div>Shim-pack Scepter C8</div>	<table><tr><td>Bonded Phase</td><td>Octyl</td></tr><tr><td>Carbon Loading</td><td>17%</td></tr><tr><td>Surface Area</td><td>360m²/g</td></tr><tr><td>Pore Size</td><td>120Å</td></tr><tr><td>Particle Size</td><td>1.9, 3, 5µm</td></tr><tr><td>pH Range</td><td>1~12</td></tr><tr><td>End-cap</td><td>Yes</td></tr><tr><td>USP Classification</td><td>L7</td></tr></table>	Bonded Phase	Octyl	Carbon Loading	17%	Surface Area	360m ² /g	Pore Size	120Å	Particle Size	1.9, 3, 5µm	pH Range	1~12	End-cap	Yes	USP Classification	L7	<div>General C8 Column:</div> <ul style="list-style-type: none">• Rapid separation for weakly polar compounds• Higher bonded density than C18 column, more effective in separation of structural isomers and analogs	<div>Pharmaceutical Analysis</div> <div>Bioanalysis</div> <div>Food / Cosmetic / Environmental Analysis</div>	<ul style="list-style-type: none">• Regular analysis with C8 column• Impurities in drug intermediates• Natural products or TCM, especially those containing alkaloids• Generic drugs bioequivalence studies• Drugs and metabolites• TCM metabolites• Regular analysis with C8 column
Bonded Phase	Octyl																			
Carbon Loading	17%																			
Surface Area	360m ² /g																			
Pore Size	120Å																			
Particle Size	1.9, 3, 5µm																			
pH Range	1~12																			
End-cap	Yes																			
USP Classification	L7																			
 <div>Shim-pack Scepter Phenyl</div>	<table><tr><td>Bonded Phase</td><td>Phenylbutyl</td></tr><tr><td>Carbon Loading</td><td>17%</td></tr><tr><td>Surface Area</td><td>360m²/g</td></tr><tr><td>Pore Size</td><td>120Å</td></tr><tr><td>Particle Size</td><td>1.9, 3, 5µm</td></tr><tr><td>pH Range</td><td>1~10</td></tr><tr><td>End-cap</td><td>Yes</td></tr><tr><td>USP Classification</td><td>L11</td></tr></table>	Bonded Phase	Phenylbutyl	Carbon Loading	17%	Surface Area	360m ² /g	Pore Size	120Å	Particle Size	1.9, 3, 5µm	pH Range	1~10	End-cap	Yes	USP Classification	L11	<div>Phenyl Column:</div> <ul style="list-style-type: none">• Powerful hydrophobic and π-π interactions• Unique selectivity for aromatic and conjugated compounds	<div>Pharmaceutical Analysis</div> <div>Bioanalysis</div> <div>Food / Cosmetic / Environmental Analysis</div>	<ul style="list-style-type: none">• Drugs and intermediates with benzene or heterocyclic rings• Drugs and metabolites• TCM metabolites• Plasticizer• Polymer additives (e.g. antioxidants)• Carcinogenic aromatic amines
Bonded Phase	Phenylbutyl																			
Carbon Loading	17%																			
Surface Area	360m ² /g																			
Pore Size	120Å																			
Particle Size	1.9, 3, 5µm																			
pH Range	1~10																			
End-cap	Yes																			
USP Classification	L11																			
 <div>Shim-pack Scepter PFPP</div>	<table><tr><td>Bonded Phase</td><td>Pentafluorophenylpropyl</td></tr><tr><td>Carbon Loading</td><td>15%</td></tr><tr><td>Surface Area</td><td>360m²/g</td></tr><tr><td>Pore Size</td><td>120Å</td></tr><tr><td>Particle Size</td><td>1.9, 3, 5µm</td></tr><tr><td>pH Range</td><td>1~8</td></tr><tr><td>End-cap</td><td>No</td></tr><tr><td>USP Classification</td><td>L43</td></tr></table>	Bonded Phase	Pentafluorophenylpropyl	Carbon Loading	15%	Surface Area	360m ² /g	Pore Size	120Å	Particle Size	1.9, 3, 5µm	pH Range	1~8	End-cap	No	USP Classification	L43	<div>Pentafluoropenylpropyl (PFPP) Column:</div> <ul style="list-style-type: none">• Powerful hydrophobic, π-π, dipole-dipole and hydrogen bonding interactions• Effective in separation of polar compounds (e.g., aromatics and halogen-containing compounds)• Trifunctional bonded stationary phases enable high reproducibility and low column bleed• High durability and stability under 100% aqueous condition	<div>Pharmaceutical Analysis</div> <div>Bioanalysis</div> <div>Food / Cosmetic / Environmental Analysis</div>	<ul style="list-style-type: none">• Highly polar compounds• Genotoxic impurities• Drugs and metabolites• TCM metabolites• Water-soluble vitamins
Bonded Phase	Pentafluorophenylpropyl																			
Carbon Loading	15%																			
Surface Area	360m ² /g																			
Pore Size	120Å																			
Particle Size	1.9, 3, 5µm																			
pH Range	1~8																			
End-cap	No																			
USP Classification	L43																			
 <div>Shim-pack Scepter Diol-HILIC</div>	<table><tr><td>Bonded Phase</td><td>Diol-HILIC</td></tr><tr><td>Carbon Loading</td><td>12%</td></tr><tr><td>Surface Area</td><td>360n²/g</td></tr><tr><td>Pore Size</td><td>120Å</td></tr><tr><td>Particle Size</td><td>1.9, 3, 5µm</td></tr><tr><td>pH Range</td><td>2~10</td></tr><tr><td>End-cap</td><td>No</td></tr><tr><td>USP Classification</td><td>L20</td></tr></table>	Bonded Phase	Diol-HILIC	Carbon Loading	12%	Surface Area	360n ² /g	Pore Size	120Å	Particle Size	1.9, 3, 5µm	pH Range	2~10	End-cap	No	USP Classification	L20	<div>High Density C18 Column:</div> <ul style="list-style-type: none">• Superior stereoselectivity, effective in separation of structural isomers and analogs• Ultra-inert, use of adsorptive ion-pairing reagent (e.g. trifluoroacetic acid, triethylamine) should be avoided in mobile phase	<div>Pharmaceutical Analysis</div> <div>Bioanalysis</div> <div>Food / Cosmetic / Environmental Analysis</div>	<ul style="list-style-type: none">• Drug products, API, intermediates and impurities• Natural products or TCM, especially those containing alkaloids• Structural isomers and analogs• Food preservatives, food dyes• Cosmetic products (preservatives, steroids, etc)
Bonded Phase	Diol-HILIC																			
Carbon Loading	12%																			
Surface Area	360n ² /g																			
Pore Size	120Å																			
Particle Size	1.9, 3, 5µm																			
pH Range	2~10																			
End-cap	No																			
USP Classification	L20																			
 <div>Shim-pack Scepter C4</div>	<table><tr><td>Bonded Phase</td><td>Butyl</td></tr><tr><td>Carbon Loading</td><td>N.D.</td></tr><tr><td>Surface Area</td><td>N.D.</td></tr><tr><td>Pore Size</td><td>300Å</td></tr><tr><td>Particle Size</td><td>1.9, 3, 5µm</td></tr><tr><td>pH Range</td><td>1~10</td></tr><tr><td>End-cap</td><td>Yes</td></tr><tr><td>USP Classification</td><td>L26</td></tr></table>	Bonded Phase	Butyl	Carbon Loading	N.D.	Surface Area	N.D.	Pore Size	300Å	Particle Size	1.9, 3, 5µm	pH Range	1~10	End-cap	Yes	USP Classification	L26	<div>Butyl Column:</div> <ul style="list-style-type: none">• Designed for separation of biomolecules such as proteins, peptides, antibodies and oligonucleotides• Excellent durability and stability even at high temperatures up to 90 °C• Good peak shapes with MS-friendly mobile phase	<div>Pharmaceutical Analysis</div> <div>Bioanalysis</div> <div>Food / Cosmetic / Environmental Analysis</div>	<ul style="list-style-type: none">• Proteins, peptides, antibodies, oligonucleotide• Therapeutic peptides• Biomarkers• Regular analysis with C4 column
Bonded Phase	Butyl																			
Carbon Loading	N.D.																			
Surface Area	N.D.																			
Pore Size	300Å																			
Particle Size	1.9, 3, 5µm																			
pH Range	1~10																			
End-cap	Yes																			
USP Classification	L26																			

Shim-pack Velox Series

Core shell particle technology enables fast analysis and high throughput

- Stable, consistent, narrow pore size distribution reduces eddy diffusion significantly, giving sharper peaks and higher sensitivity.
- Lower column backpressure, with a 2.7µm particle size column on HPLC system can achieve an UHPLC separation efficiency.
- Available in C18, SP-C18, Biphenyl, PFPP and HILIC for wide range of applications.



 <p>Shim-pack Velox C18</p>	Bonded Phase	Octadecyl	C18 Column: <ul style="list-style-type: none"> • Strong hydrophobic retention of the Shim-pack Velox phases • Ideal for analysis of complex samples • High resolution and high durability 	Pharmaceutical Analysis Bioanalysis Food / Cosmetic / Environmental Analysis	<ul style="list-style-type: none"> • Drug intermediates • Dissolution test • Acid and basic drugs and metabolites • TCM metabolites • Food preservatives, food dyes • Pesticides and veterinary drug residues • Aldehydes and ketones
	Carbon Loading	7% (2.7µm)			
	Surface Area	130m²/g (2.7µm)			
	Pore Size	90Å			
	Particle Size	1.8, 2.7, 5µm			
	pH Range	2~8			
	End-cap	Yes			
	USP Classification	L1			
 <p>Shim-pack Velox SP-C18</p>	Bonded Phase	Octadecyl	Sterically Protected C18 Column: <ul style="list-style-type: none"> • Sterically protected ligand provides extended low pH (1-3) stability, giving well-balanced retention profile and long life time • Non-encapped column, the presence of silanol group enhances hydrogen bonding, improving retention and resolution for strong acidic and basic compounds • Sharper peak shape and higher sensitivity 	Pharmaceutical Analysis Bioanalysis Food / Cosmetic / Environmental Analysis	<ul style="list-style-type: none"> • Organic synthetic process impurities • Impurities of polar drugs • Drugs and metabolites • TCM metabolites • Organic acids, nucleotides • Amino acids • PFAS
	Carbon Loading	7% (2.7µm)			
	Surface Area	130m²/g (2.7µm)			
	Pore Size	90Å			
	Particle Size	1.8, 2.7, 5µm			
	pH Range	1~8			
	End-cap	No			
	USP Classification	L1			
 <p>Shim-pack Velox Biphenyl</p>	Bonded Phase	Biphenyl	Biphenyl Column: <ul style="list-style-type: none"> • Enhanced retention of hydrophobic compounds • High selectivity for aromatic and polar compounds • Ideal for analysis of isomers and positional analogs • High durability and stability under 100% aqueous condition 	Pharmaceutical Analysis Bioanalysis Food / Cosmetic / Environmental Analysis	<ul style="list-style-type: none"> • Genotoxic impurities • Positional isomers • Positional isomers • Acid and basic drugs and metabolites • Biomarker • Metabolomics • Plasticizer • Polymer additives (e.g. antioxidants) • Veterinary drugs
	Carbon Loading	7% (2.7µm)			
	Surface Area	130m²/g (2.7µm)			
	Pore Size	90Å			
	Particle Size	1.8, 2.7, 5µm			
	pH Range	1.5~8			
	End-cap	Yes			
	USP Classification	L11			
 <p>Shim-pack Velox PFPP</p>	Bonded Phase	Pentafluorophenylpropyl	Pentafluorophenylpropyl (PFPP) Column: <ul style="list-style-type: none"> • Powerful hydrophobic, π-π, dipole-dipole and hydrogen bonding interactions • Ideal for analysis of aromatics and halogen-containing compounds, stronger retention for polar compounds than ODS column • High durability and stability under 100% aqueous condition 	Pharmaceutical Analysis Bioanalysis Food / Cosmetic / Environmental Analysis	<ul style="list-style-type: none"> • Highly polar drugs • Genotoxic impurities • Acid and basic drugs and metabolites • Isomers • Mycotoxins • Carcinogenic aromatic amines and phenols
	Carbon Loading	4% (2.7µm)			
	Surface Area	130m²/g (2.7µm)			
	Pore Size	90Å			
	Particle Size	1.8, 2.7, 5µm			
	pH Range	2~8			
	End-cap	No			
	USP Classification	L43			
 <p>Shim-pack Velox HILIC</p>	Bonded Phase	None	HILIC Column: <ul style="list-style-type: none"> • Unbonded stationary phase • Effective in separation of highly polar compounds • Excellent durability, stability and reproducibility 	Pharmaceutical Analysis Bioanalysis Food / Cosmetic / Environmental Analysis	<ul style="list-style-type: none"> • Highly polar drugs • Highly polar drugs and metabolites • Melamine, dicyandiamide • Choline, acetylcholine
	Carbon Loading	-			
	Surface Area	130m²/g			
	Pore Size	90Å			
	Particle Size	2.7µm			
	pH Range	2~8			
	End-cap	No			
	USP Classification	L3			

Shim-pack Bio Series

Solution for analysis of peptides, oligonucleotides and biopharmaceuticals

- Silica-based diol column, excellent for SEC (GFC) analysis.
- Ion-exchange column packed with low nonspecific binding hydrophilic polymer beads, wide usable pH range.
- High chemical stability HIC column, wide usable pH range.




Shim-pack Bio Diol	Particle	Silica	Size Exclusion Column: <ul style="list-style-type: none"> • Silica-based column with high mechanical stability • Suitable for analysis of mAb, oligonucleotides, carbohydrates • High resolution, fast analysis with 2µm and 3µm particles 	Pharmaceutical Analysis	<ul style="list-style-type: none"> • Therapeutic peptides • Proteins, mAb • Oligonucleotides
	Ligand	Dihydroxypropyl (Diol)			
	Particle Size	2, 3, 5µm			
	Pore Size	60Å, 120Å, 200Å, 250Å, 300Å			
	pH Range	5~7.5			
Shim-pack Bio IEX Q-NP / Q	Particle	Hydrophilic polymer	Ion Exchange Column: <ul style="list-style-type: none"> • Suitable for anion-exchange analysis • Non-porous particles offer high resolution and throughput • Porous particles offer excellent binding capacity with exceptional high efficiency, suitable for lab-scale purification 	Pharmaceutical Analysis	<ul style="list-style-type: none"> • Therapeutic peptides • Proteins • Oligonucleotides
	Ligand	-CH ₂ N ⁺ (CH ₃) ₃			
	Particle Size	3, 5µm			
	pH Range	2~12			
	Temp. Range	4-60°C			
Shim-pack Bio IEX SP-NP / SP	Particle	Hydrophilic polymer	Ion Exchange Column: <ul style="list-style-type: none"> • Suitable for cation-exchange analysis • Non-porous particles offer high resolution and throughput • Porous particles offer excellent binding capacity with exceptional high efficiency, suitable for lab-scale purification 	Pharmaceutical Analysis	<ul style="list-style-type: none"> • Therapeutic peptides • Proteins • Oligonucleotides
	Ligand	-(CH ₂) ₃ SO ₃ ⁻			
	Particle Size	3, 5µm			
	pH Range	2~12			
	Temp. Range	4-60°C			
Shim-pack Bio-HIC Butyl	Particle	Hydrophilic non-porous polymer	Hydrophobic interaction chromatography (HIC) Column: <ul style="list-style-type: none"> • Ideal for characteristic analysis of proteins and biopharmaceuticals including analysis of drug-to-antibody ratio (DAR) of antibody-drug conjugates (ADCs) • High separation capacity with low backpressure 	Pharmaceutical Analysis	<ul style="list-style-type: none"> • Antibody-Drug Conjugates (ADCs) • Proteins • Antibodies, mAb
	Ligand	Butyl			
	Particle Size	4µm			
	pH Range	2~12			
	Temp. Range	10-60°C			

Shim-pack VP/FC/XR/ Arata Series

Classic high purity silica based HPLC columns, high carbon loading and high resolution column

- High carbon loading, stronger retention, better resolution, higher durability.
- Uniform particle size distribution, high column-to-column reproducibility.
- Available in a wide range of different particle sizes and dimensions.






 <p>Shim-pack VP-ODS</p>	Bonded Phase	Octadecyl	<ul style="list-style-type: none"> • High-purity silica gel matrix with good reproducibility • Certificate of Compliance and Column Performance Report • Suitable for analysis of complex samples such as drug impurities, TCM and compound analysis 	<p>Pharmaceutical Analysis</p> <p>Food</p>	<ul style="list-style-type: none"> • API, drug products • Natural products, TCM, combination drugs • Mycotoxin
	Carbon Loading	20%			
	Surface Area	410m ² /g			
	Pore Size	120Å			
	Particle Size	5µm			
	pH Range	2~7.5			
	End-cap	Yes			
	USP Classification	L1			
 <p>Shim-pack FC-ODS</p>	Bonded Phase	Octadecyl	<ul style="list-style-type: none"> • Polymeric-bonded C18, high steric selectivity • Short column enables faster analysis 	<p>Pharmaceutical Analysis</p> <p>Food</p>	<ul style="list-style-type: none"> • API, drug products • Natural products, TCM • Mycotoxin • Azo dyes
	Carbon Loading	18%			
	Surface Area	315m ² /g			
	Pore Size	120Å			
	Particle Size	3µm			
	pH Range	1.5~9			
	End-cap	Yes			
	USP Classification	L1			
 <p>Shim-pack XR-ODS II</p>	Bonded Phase	Octadecyl	<ul style="list-style-type: none"> • High carbon load, high surface area, high resolution • Maximum pressure 60 MPa, suitable for UFLC analysis 	<p>Pharmaceutical Analysis</p>	<ul style="list-style-type: none"> • Fast analysis of API, intermediates and formulated drugs • Fast TCM analysis
	Carbon Loading	20%			
	Surface Area	470m ² /g			
	Pore Size	80Å			
	Particle Size	2.2µm			
	pH Range	2~7.5			
	End-cap	Yes			
	USP Classification	L1			
 <p>Shim-pack XR-ODS III</p>	Bonded Phase	Octadecyl	<ul style="list-style-type: none"> • High carbon load, high surface area, high resolution • Maximum pressure 100 MPa, suitable for UHPLC analysis 	<p>Pharmaceutical Analysis</p> <p>Food</p>	<ul style="list-style-type: none"> • Fast analysis of drug products, API and intermediates • Fast TCM analysis • Mycotoxin
	Carbon Loading	22%			
	Surface Area	500m ² /g (1.6µm)			
	Pore Size	75Å			
	Particle Size	1.6, 2.2µm			
	pH Range	2~7.5			
	End-cap	Yes			
	USP Classification	L1			
 <p>Shim-pack Arata C18</p>	Bonded Phase	Octadecyl	<ul style="list-style-type: none"> • Specifically designed to give unmatched peak shape for basic compounds, increased sensitivity • Ideal for basic compound analysis under low ionic strength mobile phase conditions, suitable for LC-MS analysis • Rapid column equilibrium and stable under 100% aqueous condition 	<p>Pharmaceutical Analysis</p> <p>Bioanalysis</p>	<ul style="list-style-type: none"> • Basic drugs and impurities • Therapeutic peptides • Basic drugs and their metabolites • Peptides
	Carbon Loading	17%			
	Surface Area	340m ² /g			
	Pore Size	120Å			
	Particle Size	2.2µm			
	pH Range	2~7.5			
	End-cap	Yes			
	USP Classification	L1			









Shim-pack GIST/GISS/GIS /GWS series

High purity porous spherical silica packed column, wide variety of selectivity scalable from UHPLC to preparative

- Ultra-high inertness, improve peak shapes of acidic and basic compounds.
- Uniform particle size, low backpressure, relieving load on HPLC instrument.
- Available in C18, C8, Phenyl, NH2 and other chemistries, effective over a wide range of applications.



 <p>Shim-pack GIST C18</p>	Bonded Phase	Octadecyl	General GIST C18 Column: <ul style="list-style-type: none"> • Superior inertness, provides symmetric peak shapes for strong ionic compounds • Wide selection of column dimension from analytical to preparative 	Pharmaceutical Analysis	<ul style="list-style-type: none"> • API, drug products (content / related compounds) • Drug dissolution test • Generic drugs bioequivalence studies • Drugs (acidic and basic) and metabolites • Food preservatives, food dyes • Pesticides and veterinary drug residues • PFAS
	Carbon Loading	14%		Bioanalysis	
	Surface Area	350m ² /g		Food & Environmental Analysis	
	Pore Size	100Å			
	Particle Size	2, 3, 5µm			
	pH Range	1~10			
	End-cap	Yes			
 <p>Shim-pack GIST C18-AQ</p>	Bonded Phase	Octadecyl	Optimized GIST C18 Column: <ul style="list-style-type: none"> • Optimization of C18 bonding group increases contact area with aqueous mobile phase, improving retention and resolution of polar compounds • Maintain high inertness and durability under 100% aqueous condition 	Pharmaceutical Analysis	<ul style="list-style-type: none"> • API, drug products (content / related compounds) (e.g. Amoxicillin) • Generic drugs bioequivalence studies • Drugs (acidic and basic) and metabolites
	Carbon Loading	13%		Bioanalysis	
	Surface Area	350m ² /g			
	Pore Size	100Å			
	Particle Size	1.9, 3, 5µm			
	pH Range	1~10			
	End-cap	Yes			
 <p>Shim-pack GISS C18</p>	Bonded Phase	Octadecyl	GISS C18 Column: <ul style="list-style-type: none"> • Low carbon load, low surface area, rapid separation with symmetrical peaks • Ideal for high speed LC/MS/MS analysis 	Pharmaceutical Analysis	<ul style="list-style-type: none"> • API, drug products (content / related compounds) • Biological samples
	Carbon Loading	9%		Bioanalysis	
	Surface Area	200m ² /g			
	Pore Size	200Å			
	Particle Size	1.9, 3, 5µm			
	pH Range	1~10			
	End-cap	Yes			
 <p>Shim-pack GIST C8</p>	Bonded Phase	Octyl	GIST C8 Column: <ul style="list-style-type: none"> • High inertness, effective in separation of acidic and basic compounds with symmetrical peaks 	Pharmaceutical Analysis	<ul style="list-style-type: none"> • API, drug products (content / related compounds) • Generic drugs bioequivalence studies • Drugs (acidic and basic) and metabolites
	Carbon Loading	8%		Bioanalysis	
	Surface Area	350m ² /g			
	Pore Size	100Å			
	Particle Size	2, 3, 5µm			
	pH Range	1~10			
	End-cap	Yes			
 <p>Shim-pack GIST Phenyl</p>	Bonded Phase	Phenyl	Phenyl Column: <ul style="list-style-type: none"> • Phenyl group directly bonded to silica, providing stronger π-π interactions • Suitable for aromatics or heterocyclic and electronegative compounds • Non-encapped column, the presence of silanol group provides hydrogen bonding, improving retention and resolution of acidic and basic compounds 	Pharmaceutical Analysis	<ul style="list-style-type: none"> • Electronegative and benzene / heterocyclic ring containing drug and drug impurities • Drugs (acidic and basic) and metabolites
	Carbon Loading	10%		Bioanalysis	
	Surface Area	350m ² /g			
	Pore Size	100Å			
	Particle Size	2, 3, 5µm			
	pH Range	2~7.5			
	End-cap	No			
 <p>Shim-pack GIST Phenyl-Hexyl</p>	Bonded Phase	Phenylhexyl	Phenyl-Hexyl Column: <ul style="list-style-type: none"> • Phenylhexyl-bonded column, provides π-π and hydrophobic interactions • Ideal for analysis of aromatics, heterocyclic and other electronegative compounds 	Pharmaceutical Analysis	<ul style="list-style-type: none"> • Electronegative and benzene / heterocyclic ring containing drug and drug impurities • Benzene or heterocyclic ring containing drugs and its metabolites • Plasticizer • PAHs • Phenols
	Carbon Loading	9%		Bioanalysis	
	Surface Area	350m ² /g		Food & Environmental Analysis	
	Pore Size	100Å			
	Particle Size	3, 5µm			
	pH Range	1~10			
	End-cap	Yes			
 <p>Shim-pack GIST NH2</p>	Bonded Phase	Aminopropyl	NH2 Column: <ul style="list-style-type: none"> • Highly versatile column bonded with aminopropyl group • Effective in separation of sugars, amino acids and polar compounds under HILIC mode • Superior stability, reduced column bleed under weakly acidic eluent 	Pharmaceutical Analysis	<ul style="list-style-type: none"> • Carbohydrates (mono- and disaccharides), amino acids • Polar compounds and metabolites • Monosaccharides and disaccharides
	Carbon Loading	7%		Bioanalysis	
	Surface Area	350m ² /g		Food & Environmental Analysis	
	Pore Size	100Å			
	Particle Size	3, 5µm			
	pH Range	2~7.5			
	End-cap	No			
 <p>Shim-pack GIS C18</p>	Bonded Phase	Octadecyl	GIS C18 Column: <ul style="list-style-type: none"> • High-purity fully porous silica matrix, high surface area, high loading capacity • Ideal for analysis of synthetic drugs and natural products 	Pharmaceutical Analysis	<ul style="list-style-type: none"> • API, drug products (content / related compounds) • Chinese herbal medicine (raw material and prescriptions) • Preparative isolation and purification of synthetic drugs and natural products
	Carbon Loading	15%			
	Surface Area	450m ² /g			
	Pore Size	100Å			
	Particle Size	2, 3, 4, 5, 10µm			
	pH Range	2~7.5			
	End-cap	Yes			
	USP Classification	L1			

 Shim-pack GWS C18	Bonded Phase	Octadecyl	GWS C18 Column: <ul style="list-style-type: none"> High-purity fully porous silica matrix, low carbon load, improved retention of polar compounds, faster elution of hydrophobic compounds Ideal for analysis of multi-component samples with wide range of polarities 	Pharmaceutical Analysis	<ul style="list-style-type: none"> API, drug products (content / related compounds)
	Carbon Loading	9.5%			
	Surface Area	450m²/g			
	Pore Size	100Å			
	Particle Size	5µm			
	pH Range	2~7.5			
	End-cap	Yes			
	USP Classification	L1			
 Shim-pack GIS C18-P	Bonded Phase	Octadecyl	Polymeric C18 Column: <ul style="list-style-type: none"> High steric selectivity for separation of planar and non-planar compounds High carbon load, high surface area, strong retention force 	Pharmaceutical Analysis Bioanalysis Food & Environmental Analysis	<ul style="list-style-type: none"> Drug products, API, intermediates and impurities Natural products or TCM, especially those containing alkaloids Structural isomers and analogs Food preservatives, food dyes, PAHs Pesticides and veterinary drug residues
	Carbon Loading	29%			
	Surface Area	450m²/g			
	Pore Size	100Å			
	Particle Size	3, 5µm			
	pH Range	2~7.5			
	End-cap	No			
	USP Classification	L1			
 Shim-pack GIS RP-Shield	Bonded Phase	Octadecyl	Polar embedded C18 Column: <ul style="list-style-type: none"> Extreme base deactivation enables superior peak shape for acids Unique selectivity through hydrogen bonding interactions Maintain high inertness and durability under 100% aqueous condition 	Pharmaceutical Analysis Bioanalysis Food & Environmental Analysis	<ul style="list-style-type: none"> API, drug products (content / related compounds) Drugs (acidic and basic) and metabolites Food preservatives, food dyes Organic acids
	Carbon Loading	9%			
	Surface Area	450m²/g			
	Pore Size	100Å			
	Particle Size	5µm			
	pH Range	2~7.5			
	End-cap	No			
	USP Classification	L1			
 Shim-pack GIST PFPP	Bonded Phase	Pentafluorophenylpropyl	Pentafluoropropyl (PFPP) Column: <ul style="list-style-type: none"> End-capped column, better selectivity than non-endcapped column Provides hydrophobic, π-π, dipole-dipole and hydrogen bonding interactions Ideal for analysis of polar compounds, aromatics and halogen-containing compounds 	Pharmaceutical Analysis Bioanalysis Food & Environmental Analysis	<ul style="list-style-type: none"> Highly polar drugs Genotoxic impurities Drugs and metabolites analysis Isomers Mycotoxins Carcinogenic aromatic amines and phenols
	Carbon Loading	10%			
	Surface Area	350m²/g			
	Pore Size	100Å			
	Particle Size	3, 5µm			
	pH Range	2~7.5			
	End-cap	Yes			
	USP Classification	L43			
 Shim-pack GIST Amide	Bonded Phase	Amide (carbamoyl)	Amide-type HILIC Column: <ul style="list-style-type: none"> Stronger retention for polar compounds High stability and durability, short column equilibration time 	Pharmaceutical Analysis Bioanalysis Food & Environmental Analysis	<ul style="list-style-type: none"> Polar compounds Carbohydrates, amino acids Polar compounds and metabolites Glycans Melamine, cyanuric acid Amino acids, water-soluble vitamins
	Carbon Loading	15%			
	Surface Area	350m²/g			
	Pore Size	100Å			
	Particle Size	3, 5µm			
	pH Range	2~8.5			
	End-cap	No			
	USP Classification	L68			
 Shim-pack GIS CN	Bonded Phase	Cyanopropyl	Cyanopropyl Column: <ul style="list-style-type: none"> Provides π-π and hydrogen bonding interactions Compatible with both reversed phase and normal phase conditions Applicable in 100% aqueous condition 	Pharmaceutical Analysis Bioanalysis Food & Environmental Analysis	<ul style="list-style-type: none"> Highly polar drugs Highly polar drugs and metabolites Amino acids, organic acids Water-soluble vitamins
	Carbon Loading	14%			
	Surface Area	450m²/g			
	Pore Size	100Å			
	Particle Size	3, 5µm			
	pH Range	2~7.5			
	End-cap	No			
	USP Classification	L10			
 Shim-pack GIS HILIC	Bonded Phase	Diol	HILIC Column: <ul style="list-style-type: none"> Dihydroxypropyl (Diol) bonded stationary phase Excellent peak shape for basic and neutral polar compounds 	Pharmaceutical Analysis Bioanalysis Food / Cosmetic / Environmental Analysis	<ul style="list-style-type: none"> Highly polar drugs Highly polar drugs and metabolites Allantoin, melamine, cyanoguanidine Taurine, inositol
	Carbon Loading	20%			
	Surface Area	450m²/g			
	Pore Size	100Å			
	Particle Size	3, 5µm			
	pH Range	2~7.5			
	End-cap	No			
	USP Classification	L20			
 GIS SIL	Bonded Phase	None	Silica Column: <ul style="list-style-type: none"> High quality pure silica gel compatible with normal phase condition Strongly retain basic compounds due to the silanol groups of the silica gel 	Pharmaceutical Analysis Bioanalysis Food / Cosmetic / Environmental Analysis	<ul style="list-style-type: none"> API, drug products (content / related compounds) Drug dissolution test Steroid hormones Biofuel, FAME Fat-soluble vitamins Cosmetic preservatives
	Carbon Loading	-			
	Surface Area	450m²/g			
	Pore Size	100Å			
	Particle Size	3, 5µm			
	pH Range	2~7.5			
	End-cap	No			
	USP Classification	L3			

Shim-pack NovaCore

Core shell organosilica hybrid particle technology enables fast analysis with wide pH range

- Stable, consistent, narrow pore size distribution reduces eddy diffusion significantly, giving sharper peaks and higher sensitivity.
- Lower column backpressure, with a 2.6µm particle size column on HPLC system can achieve an UHPLC separation efficiency.
- Organosilica hybrid base material offer wide pH tolerance range (pH 1 to 12).



Shim-pack
NovaCore
C18-HB

Bonded Phase	Octadecyl
Carbon Loading	11%
Surface Area	200m ² /g
Pore Size	100Å
Particle Size	1.7, 2.6, 5µm
pH Range	1~12
End-cap	Yes
USP Classification	L1

C18 Column:

- Ideal for method development over wide pH range
- High resolution and high durability

Pharmaceutical Analysis

- Drug intermediates
- Acid and basic drugs and metabolites

Food / Cosmetic / Environmental Analysis

- Food preservatives, food dyes
- Synthetic dye
- Cosmetic products (preservatives, steroids, etc.)
- Pesticides and veterinary drug residues

Mobile Phase Cleaner

Mobile phase cleaner for UHPLC/HPLC


- Consistently traps impurities in mobile phases, effectively removes ghost peaks.
- Ultra high surface area activated charcoal beads, large adsorption capacity and long lifetime.



Ghost Trap DS

Ghost Trap DS Handling and Precautions

- Ghost Trap DS is installed between the gradient mixer and the autosampler.
- Cleans up impurities from piping and gradient mixer in addition to those from the mobile phase.
- Not recommended for analysis with ion-pairing reagent (e.g. triethylamine, trifluoroacetic acid, sodium alkyl sulfonate).
- Suitable for LC-UV analysis, not recommend to use with MS detector.
- Thoroughly rinse the system and Ghost Trap DS with 100% acetonitrile before loading the column.
- Product lifetime varies according to LC analysis conditions and mobile phase used.
- Note that some impurities may not be removed.



GLC Suction Filter 3

HPLC Suction filter with mobile phase cleaning function

- Removes contaminants from the mobile phase and prevents airborne contaminants from entering the mobile phase.
- Does not affect gradient delay to install this product in front of the pump.
- Can be used immediately by replacing the existing suction filter, as conditioning before first use is generally unnecessary.




Fig. 1. Structure of the GLC Suction Filter 3

GLC Suction Filter 3 Handling and Precautions

- These products do not guarantee compatibility with all instruments. Please confirm the outer diameter of the suction tube before use.
- When using ion-pair reagents in the mobile phase, the ion-pair reagents may be adsorbed to these products and may affect the retention time and peak shape.
- Bleeding may occur during the use of mass spectrometers, depending on the conditions of the mobile phase. In addition, bleeding originating from this product may occur during initial use. In such cases, please condition the product.
- These products do not remove all contaminants.
- When the activated carbon reaches breakthrough, "ghost peaks" derived from mobile phase contaminants may be detected. In this case, replace the product.
- The recommended replacement period is after passing 60 L of mobile phase or 1 year after use.




Fig. 2. Attachment position

Column Selection Guide for Different LC Platforms

- Column configuration should be optimise according to the LC system used. Smaller column particle size, or reduced column volume (ID and/or length of column) increases the necessity for a lower dispersion system.
- Choosing the optimal column configuration for your LC system allows you to achieve improved chromatography. The following table summarizes the starting recommendations of column configuration for each LC system.



Nexera Lite



LC-2070



LC-2080




Nexera XS



Nexera X4

LC System	HPLC	UHPLC-like	UHPLC
Particle Size	2.7, 3, 5µm	2.7, 3µm	1.6, 1.8, 1.9, 2, 2.2µm
Column I.D.	4.6mm (3.0mm)	3.0mm (2.1mm)	2.1mm
Column Length	100 - 250mm	50 - 100mm	≤ 150mm




Method Development Column Selection Flowchart

First Choice Columns

- Scepter C18-120 (Operating pH range: 1 to 12)
- GIST C18 (Operating pH range: 1 to 10)

Red: Fully porous silica columns
Green: Fully porous organic silica hybrid columns
Blue: Surface porous silica columns (Core-shell silica columns)

When satisfactory analysis is not possible even if the analysis conditions (including mobile phase solvent, solvent ratio, mobile phase pH conditions, and temperature) are changed



*1 This is not a description of all of the Shim-pack columns.

*2 The separation selectivity differs significantly even for columns in the same group. Actual column selection will be more complicated than shown in this chart depending on the analytes and conditions.

The analysis results in this column selection chart are not guaranteed, so use this only as a reference.



Shimadzu (Asia Pacific) Pte Ltd.

www.shimadzu.com.sg/an/

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.